



COMPLETE TESTING SOLUTIONS

Where a negative is a positive.

(480) 568-8222

service@ctsdrugtesting.com

Manufactured by: HEALGEN

# UA Cups

Our Screens are CLIA Waived,  
FDA 510(k) Cleared, and adhere to  
SAMSHA cutoff levels.



## 5 Panel

Item #	Configuration	Available with Adulterants/A3 (OX/SG/pH)
154	THC, COC, OPI2000, AMP, PCP	X
254	THC, COC, OPI2000, AMP, MET	X
654	THC, COC, OPI2000, MET, BZO	X
855	COC, MOP300, AMP, MET, OXY	
955	COC, MOP300, BZO, BUP, MTD	

## 6 Panel

164	THC, COC, OPI2000, AMP, MET, PCP	
264	THC, COC, OPI2000, AMP, MET, BZO	X
265	THC, COC, MOP300, AMP, MET, BZO	
564	THC, COC, OPI2000, MET, BZO, OXY	

## 7 Panel

274	THC, COC, OPI2000, AMP, MET, BZO, OXY	X
474	THC, COC, OPI2000, AMP, MET, PCP, MDMA	
874	THC, COC, OPI2000, AMP, PCP, BAR, MDMA	
975	COC, MOP300, AMP, BZO, BAR, OXY, MTD	

# 10 Panel

Item #	Configuration	Available with Adulterants/A3 (OX/SG/pH)
<b>1104k</b>	K2, THC, COC, OPI2000, AMP, MET, BZO, BUP, OXY, MDMA	
<b>2104</b>	THC, COC, OPI2000, AMP, MET, BZO, BAR, OXY, MDMA, PPX	
<b>3104</b>	THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, MDMA, MTD	<b>X</b>
<b>4104</b>	THC, COC, OPI2000, AMP, MET PCP, BZO, BAR, OXY, MTD	<b>X</b>
<b>7104</b>	THC, COC, OPI2000, MET, BZO, BAR, BUP, OXY, MTD, TCA	<b>X</b>
<b>7105</b>	THC, COC, MOP300, MET, BZO, BAR, BUP, OXY, MTD, TCA	
<b>8104</b>	THC, COC, OPI2000, AMP, MET, BZO, BUP, OXY, MDMA, MTD	

# 12 Panel

<b>2125</b>	THC, COC, MOP300, AMP, MET, BZO, BAR, OXY, MDMA, MTD, TCA, PPX	
<b>2125k</b>	K2, THC, COC, MOP300, AMP, MET, BZO, BAR, BUP, OXY, MDMA, MTD	
<b>3124</b>	THC, COC, OPI2000, AMP, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA	<b>X</b>
<b>3125</b>	THC, COC, MOP300, AMP, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA	
<b>4124</b>	THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, BUP, OXY, MTD, PPX	<b>X</b>
<b>4125</b>	THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MTD, PPX	<b>X</b>
<b>5125</b>	THC, MET, COC, MOP, BZO, AMP, OXY, MDMA, FEN, MTD, BAR, BUP	<b>X</b>
<b>6124</b>	THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD	<b>X</b>
<b>6125</b>	THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD	<b>X</b>
<b>9125</b>	THC, COC, MOP300, AMP, MET, BZO, BAR, BUP, OXY, MDMA, MTD, TCA	

# 13 Panel



Item #	Configuration	Available with Adulterants/A3 (OX/SG/pH)
2135	THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA	X
3135SETF	FEN, TRA, ETG300, THC, COC, MOP300, AMP, MET, BZO, BUP, OXY, MDMA, MTD	
4135SEFKT	K2, FEN, EtG, TRA, THC, COC150, OPI300, AMP300, BZO, BUP, OXY, MDMA500, EDDP	
6135	ETG, THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD	X
7135EFT	AMP, BZO, ETG, COC, MDMA, MTD, MAMP, FEN, OPI300, OXY, THC, BUP, TRA	

# 14 Panel



1144	THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA, PPX	X
1145	THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA, PPX	X
2145	THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA, EDDP	X
6145EF	FEN, ETG, THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD	
8145EF	AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FEN	

# 15 Panel



2155EF	FEN, ETG, THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA	
6155E3FT	FEN, TRA, ETG300, THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD	

# 16 Panel



6155E3FT	FEN, TRA, ETG300, THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD	
6164EFKT	K2, FEN, TRA, ETG, THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD	

# EZ SPLIT KEY CUPS



5-Panel	<b>DOA-1257-019</b>	AMP / COC / mAMP / OPI / THC
5-Panel	<b>DOA-3257-019</b>	COC / mAMP / OPI / PCP / THC
6-Panel	<b>DOA-167-321-01</b>	COC / mAMP / MDMA / OPI / OXY / THC
6-Panel	<b>DOA-2167-019</b>	AMP / BZO / COC / mAMP / OPI / THC
12-Panel	<b>DOA-1127-041-019</b>	AMP / BAR / BUP / BZO / COC / mAMP / MDMA / MTD / OPI / OXY / PPX / THC

- Locking key mechanism to prevent tampering & adulteration.
- Easily Photocopy Results with flat-side of cup.
- Ready for lab send-off.



Need assistance or volume-based pricing?

Setting up an account with CTS is easy, just give us a quick call.

**(480) 568-8222**

email: [service@ctsdrugtesting.com](mailto:service@ctsdrugtesting.com)





# Cut-off Levels

Test	Cut-off (ng/mL)	Detection Time	Certifications
Amphetamine (AMP)	300, 500 or 1000	2 hours to 4 days	FDA CLIA Waived (500, 1000)
Barbituates (BAR)	200 or 300	2 hours to 3 days	FDA CLIA Waived (300)
Benzodiazepines (BZO)	200, 300 and 500	2 hours to 4 days	FDA CLIA Waived (300)
Buprenorphine (BUP)	5 or 10	2 hours to 3 days	FDA CLIA Waived (10)
Cocaine (COC)	100, 150, 200 or 300	1 hours to 4 days	FDA CLIA Waived (150, 300)
Cannabinoids (THC)	25, 30, 50, 150 or 200	2 hours to 40 days	FDA CLIA Waived (50)
Methadone (MTD)	200 or 300	3 hours to 3 days	FDA CLIA Waived (300)
Methamphetamine (MET)	300, 500 or 1000	2 hours to 4 days	FDA CLIA Waived (500, 1000)
Morphine (MOP)	100, 200 or 300	2 hours to 3 days	FDA CLIA Waived (300)
Opiates (OPI)	2000	2 hours to 3 days	FDA CLIA Waived
Oxycodone (OXY)	100	1 hours to 2 days	FDA CLIA Waived
Phencyclidine (PCP)	25	4 hours to 14 days	FDA CLIA Waived
D,L-Methylenedioxy-Methamphetamine(MDMA)	300, 500 or 1000	2 hours to 4 days	FDA CLIA Waived (500)
Methadone Metabolite (EDDP)	100 or 300	1 hours to 3 days	FDA CLIA Waived (100)
Nortriptyline (TCA)	1000	8 hours to 7 days	FDA CLIA Waived
Propoxyphene (PPX)	300	1 hours to 3 days	FDA CLIA Waived

7-Aminoclonazepam (7-ACL)	Available Soon	Available Soon	Forensic Use Only
Alcohol (ALC)	0.02% BAC	N/A	Forensic Use Only
Caffiene	Available Soon	Available Soon	Forensic Use Only
Cotinine(COT)	100 or 200	2 hours to 3 days	Forensic Use Only
Ethyl Glucuronide (ETG)	300, 500 or 1000	1 hours to 3 days	Forensic Use Only
Fentanyl (FEN)	200 or 300	N/A	Forensic Use Only
Synthetic Marijuana (K2)	25 or 50	1 hour to 2 days	Forensic Use Only
Ketamine (KET)	100, 300, 500 or 1000	Up to 2 days	Forensic Use Only
LSD	20	1 to 2 days	Forensic Use Only
Methcathinone (MCAT500)	500	N/A	Forensic Use Only
Methylphenidate (MPD)	300	1 hour to 3 days	Forensic Use Only
Methqualone (MQL)	300	2 hours to 4 days	Forensic Use Only
Bath Salts (MDPV)	1000	1 to 2 days	Forensic Use Only
Tramadol (TRA)	100, 200	1 hour to 2 days	Forensic Use Only
6-monoacetylmorphine (6-MAM)	10	24 hours	Forensic Use Only

# EZ Split Key Cup Cut-off Levels

Drug Name	ABBR	Calibrator	Standard Cutoff Levels Urine	Standard Cutoff Levels Saliva	Approx. Detection Time Urine
Alcohol	ALC	Ethanol	>0.04%	>0.02%	1-12 hours
Amphetamine	AMP	d-Amphetamine	500 & 1000 ng/ml	50 ng/ml	2-4 Days
Barbiturates	BAR	Secobarbital	300 ng/ml	50 ng/ml	1-3 Days
Benzodiazepines	BZO	Oxazepam	300 ng/ml	10 ng/ml	3-5 Days
Buprenorphine	BUP	Buprenorphine	10 ng/ml	5 ng/ml	3-5 Days
Cocaine	COC	Benzoylcegonine	150 & 300 ng/ml	20 ng/ml	2-4 Days (Single) > 2 Weeks (Long Term)
Cotinine	COT	Cotinine	200 ng/ml	30 ng/ml	1-3 Days
Ecstasy	MDMA	D,L-Methylenedioxy-methamphetamine	500 ng/ml	NA	1-4 Days

**Please Note: Cut-off charts give approximate detection periods for each substance. The ranges depend on amount and frequency of use, metabolic rate, body mass, age, overall health, drug tolerance, urine pH, and other additional factors.**

<b>Morphine</b>	MOP	Morphine	300 ng/ml	40 ng/ml	1-3 Days
<b>Opiates</b>	OPI	Morphine	2000 ng/ml		
<b>Oxycodone</b>	OXY	Oxycodone	100 ng/ml	20 ng/ml	1-3 Days
<b>Phencyclidine</b>	PCP	Phencyclidine	25 ng/ml	10 ng/ml	2-13 Days
<b>Tramadol</b>	TRA	Tramadol	200 ng/ml	NA	1-3 Days
<b>Tricyclic Antidepressants</b>	TCA	Nortriptyline	1000 ng/ml	NA	2-5 Days
<b>Propoxyphene</b>	PPX	Propoxyphene	300 ng/ml	NA	1-7 Days
<b>Ethyl Glucuronide</b>	ETG	Ethyl Glucuronide	300 & 500 ng/ml	NA	Up to 80 Hours
<b>Fentanyl</b>	FEN	Norfentanyl	200 & 300 ng/ml	10 ng/ml	3 Days
<b>Gabapentin</b>	GABA	Gabapentin	1000 ng/ml	N/A	1-2 Days
<b>K2 / Spice</b>	K2/Spice	JWH-073/JWH-018	50 ng/ml	5 ng/ml	1-3 Days
<b>Ketamine</b>	KET	Ketamine	1000 ng/ml	N/A	1-3 Days
<b>Marijuana</b>	THC	11-nor- $\Delta^9$ -THC-9 COOH	50 ng/ml	12 ng/ml	2-5 Days (Single) 14-21 days (Long Term)
<b>Methadone</b>	MTD	Methadone	50 ng/ml	30 ng/ml	1-7 days
<b>Methamphetamine</b>	MAMP	D-Methamphetamine	500 & 1000 ng/ml	50 ng/ml	2-4 Days