



COMPLETE TESTING SOLUTIONS

Where a negative is a positive.

(480) 568-8222

service@ctsdrugtesting.com

Manufactured by: HEALGEN

# Dip Cards

Our Screens are CLIA Waived,  
FDA 510(k) Cleared, and adhere to  
SAMSHA cutoff levels.



# Single Panels



|  |
|--|
| <b>AMP/Amphetamines</b>                    |
| <b>BAR/Barbiturates</b>                    |
| <b>BUP/Buprenorphine</b>                   |
| <b>BZO/Benzodiazepine</b>                  |
| <b>CLO/Clonazepam (Klonopin) 100 ng/mL</b> |
| <b>COC/Cocaine</b>                         |
| <b>COT/Cotinine (Nicotine Metabolite)</b>  |
| <b>ETG500 (Alcohol Metabolite)</b>         |
| <b>ETG300 (Alcohol Metabolite)</b>         |
| <b>FTN/Fentanyl</b>                        |
| <b>K2 (Synthetic Marijuana) 50 ng/mL</b>   |
| <b>K2 (Synthetic Marijuana) 25 ng/mL</b>   |
| <b>KET/Ketamine 1000 ng/mL</b>             |
| <b>KRA/Kratom 250 ng/ml</b>                |
| <b>LSD Dip Card</b>                        |
| <b>MET/Methamphetamine</b>                 |

# Single Panels

|                               |
|-------------------------------|
| MDMA/Ecstasy                  |
| MOP300/Morphine               |
| MPD/Methylphenidate           |
| MQL (Methaqualone) 300 ng/mL  |
| MTD/Methadone                 |
| MDPV (Bath Salts)             |
| OPI2000/Opiate                |
| OXY/Oxycodone                 |
| PCP/Phencyclidine             |
| PGB/Pregabalin 1000 ng/mL     |
| PPX/Propoxyphene              |
| TCA/Tricyclic Antidepressants |
| THC/Marijuana                 |
| TRA/Tramadol                  |
| 6-MAM (6-Monoacetylmorphine)  |

## 3 Panel

| Item # | Configuration     | Available with Adulterants/A3 (OX/SG/pH) |
|--------|-------------------|--|
| 134    | THC, COC, OPI2000 |  |
| 234    | THC, COC, MET     |  |

# 5 Panel

| Item # | Configuration               | Available with Adulterants/A3 (OX/SG/pH) |
|--------|-----------------------------|--|
| 154    | THC, COC, OPI2000, AMP, PCP | X  |
| 254    | THC, COC, OPI2000, AMP, MET | X  |
| 255    | THC, COC, MOP300, AMP, MET  |  |
| 354    | THC, COC, OPI2000, MET, OXY |  |
| 554    | THC, COC OPI2000, MET, MDMA |  |
| 654    | THC, COC, OPI2000, MET, BZO | X  |
| 754    | THC, COC, OPI2000, AMP, BZO |  |

# 6 Panel

|     |                                  |  |
|-----|----------------------------------|--|
| 164 | THC, COC, OPI2000, AMP, MET, PCP |  |
| 264 | THC, COC, OPI2000, AMP, MET, BZO |  |
| 265 | THC, COC, MOP300, AMP, MET, BZO  |  |
| 564 | THC, COC, OPI2000, MET, BZO, OXY |  |
| 565 | THC, COC, MOP300, MET, BZO, OXY  |  |
| 664 | THC, COC, OPI2000, AMP, MET, OXY |  |
| 865 | COC, MOP300, BZO, BUP, OXY, MTD  |  |

# 7 Panel

|     |                                       |  |
|-----|---------------------------------------|--|
| 274 | THC, COC, OPI2000, AMP, MET, BZO, OXY |  |
| 975 | COC, MOP300, AMP, BZO, BAR, OXY, MTD  |  |

# 8 Panel

|       |  |  |
|-------|--|--|
| 184   | THC, COC, OPI2000, MET, AMP, PCP, BUP, BZO   |  |
| 185ET | ETG500, THC, COC, MOP300, MET, BZO, BUP, OXY |  |
| 285   | THC, COC, MOP300, AMP, MET, BZO, BUP, OXY    |  |

# 9 Panel



| Item # | Configuration                                   | Available with Adulterants/A3 (OX/SG/pH) |
|--------|---|--|
| 194    | THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, MTD |  |
| 295    | THC, COC, MOP300, AMP, BZO, BUP, OXY, MTD, PPX  |  |

# 10 Panel



|           |   |   |
|-----------|---|---|
| 1104      | THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, MTD, TCA  | X |
| 1105      | THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, MTD, TCA   |   |
| 2105      | THC, COC, MOP300, AMP, MET, BZO, BAR, OXY, MTD, TCA   |   |
| 3104      | THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, MDMA, MTD |   |
| 4104      | THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, OXY, MTD  | X |
| 4105      | THC, COC, MOP300, AMP, MET, PCP, OXY, BZO, BAR, MTD   |   |
| 6104      | THC, COC, OPI2000, AMP, BZO, BAR, BUP, OXY, MTD, TCA  |   |
| 7104      | THC, COC, OPI2000, MET, BZO, BAR, BUP, OXY, MTD, TCA  | X |
| 7105      | THC, COC, MOP300, MET, BZO, BAR, BUP, OXY, MTD, TCA   |   |
| 8104      | THC, COC, OPI2000, AMP, MET, BZO, BUP, OXY, MDMA, MTD |   |
| 16105E3GT | GAB, ETG300, TRA, THC, COC, MOP, MET, BZO, BUP, OXY   |   |

# 12 Panel

| Item # | Configuration   | Available with Adulterants/A3 (OX/SG/pH) |
|--------|---|--|
| 2125   | THC, COC, MOP300, AMP, MET, BZO, BAR, OXY, MDMA, MTD, TCA, PPX  |  |
| 3124   | THC, COC, OPI2000, AMP, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA |  |
| 5125   | THC, COC, MOP300, AMP, MET, BAR, BZO, MDMA, MTD, OXY, PCP, TCA  |  |
| 5125   | THC, MET, COC, MOP, BZO, AMP, OXY, MDMA, FEN, MTD, BAR, BUP     |  |
| 6124   | THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD |  |
| 6125   | THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD  | X  |

# 14 Panel

|      |   |  |
|------|---|--|
| 1144 | THC, COC, OPI2000, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA, PPX |  |
| 1145 | THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA, PPX  |  |
| 2145 | THC, COC, MOP300, AMP, MET, PCP, BZO, BAR, BUP, OXY, MDMA, MTD, TCA, EDDP |  |

Need assistance or volume-based pricing?

Setting up an account with CTS is easy, just give us a quick call.

**(480) 568-8222**

email: [service@ctsdrugtesting.com](mailto:service@ctsdrugtesting.com)





# Cut-off Levels

| Test                                     | Cut-off (ng/mL)        | Detection Time     | Certifications              |
|--|------------------------|--------------------|-----------------------------|
| Amphetamine (AMP)                        | 300, 500 or 1000       | 2 hours to 4 days  | FDA CLIA Waived (500, 1000) |
| Barbituates (BAR)                        | 200 or 300             | 2 hours to 3 days  | FDA CLIA Waived (300)       |
| Benzodiazepines (BZO)                    | 200, 300 and 500       | 2 hours to 4 days  | FDA CLIA Waived (300)       |
| Buprenorphine (BUP)                      | 5 or 10                | 2 hours to 3 days  | FDA CLIA Waived (10)        |
| Cocaine (COC)                            | 100, 150, 200 or 300   | 1 hours to 4 days  | FDA CLIA Waived (150, 300)  |
| Cannabinoids (THC)                       | 25, 30, 50, 150 or 200 | 2 hours to 40 days | FDA CLIA Waived (50)        |
| Methadone (MTD)                          | 200 or 300             | 3 hours to 3 days  | FDA CLIA Waived (300)       |
| Methamphetamine (MET)                    | 300, 500 or 1000       | 2 hours to 4 days  | FDA CLIA Waived (500, 1000) |
| Morphine (MOP)                           | 100, 200 or 300        | 2 hours to 3 days  | FDA CLIA Waived (300)       |
| Opiates (OPI)                            | 2000                   | 2 hours to 3 days  | FDA CLIA Waived             |
| Oxycodone (OXY)                          | 100                    | 1 hours to 2 days  | FDA CLIA Waived             |
| Phencyclidine (PCP)                      | 25                     | 4 hours to 14 days | FDA CLIA Waived             |
| D,L-Methylenedioxy-Methamphetamine(MDMA) | 300, 500 or 1000       | 2 hours to 4 days  | FDA CLIA Waived (500)       |
| Methadone Metabolite (EDDP)              | 100 or 300             | 1 hours to 3 days  | FDA CLIA Waived (100)       |
| Nortriptyline (TCA)                      | 1000                   | 8 hours to 7 days  | FDA CLIA Waived             |
| Propoxyphene (PPX)                       | 300                    | 1 hours to 3 days  | FDA CLIA Waived             |

|                              |                       |                   |                   |
|------------------------------|-----------------------|-------------------|-------------------|
| 7-Aminoclonazepam (7-ACL)    | Available Soon        | Available Soon    | Forensic Use Only |
| Alcohol (ALC)                | 0.02% BAC             | N/A               | Forensic Use Only |
| Caffeine                     | Available Soon        | Available Soon    | Forensic Use Only |
| Cotinine(COT)                | 100 or 200            | 2 hours to 3 days | Forensic Use Only |
| Ethyl Glucuronide (ETG)      | 300, 500 or 1000      | 1 hours to 3 days | Forensic Use Only |
| Fentanyl (FEN)               | 200 or 300            | N/A               | Forensic Use Only |
| Synthetic Marijuana (K2)     | 25 or 50              | 1 hour to 2 days  | Forensic Use Only |
| Ketamine (KET)               | 100, 300, 500 or 1000 | Up to 2 days      | Forensic Use Only |
| LSD                          | 20                    | 1 to 2 days       | Forensic Use Only |
| Methcathinone (MCAT500)      | 500                   | N/A               | Forensic Use Only |
| Methylphenidate (MPD)        | 300                   | 1 hour to 3 days  | Forensic Use Only |
| Methqualone (MQL)            | 300                   | 2 hours to 4 days | Forensic Use Only |
| Bath Salts (MDPV)            | 1000                  | 1 to 2 days       | Forensic Use Only |
| Tramadol (TRA)               | 100, 200              | 1 hour to 2 days  | Forensic Use Only |
| 6-monoacetylmorphine (6-MAM) | 10                    | 24 hours          | Forensic Use Only |

**Please Note: Cut-off charts give approximate detection periods for each substance. The ranges depend on amount and frequency of use, metabolic rate, body mass, age, overall health, drug tolerance, urine pH, and other additional factors.**