



COMPLETE TESTING SOLUTIONS

Where a negative is a positive.

(480) 568-8222

service@ctsdrugtesting.com

Our Screens are CLIA Waived,
FDA 510(k) Cleared, and adhere to
SAMSHA cutoff levels.

Oral Fluid Saliva Swabs



Device E with Saliva Saturation Indicator



Device F with Saliva Saturation Indicator
(Optional)

5 Panel

Item #	Configuration & Cut-off	Saliva Saturation Indicator
7455ESI	THC12, COC20, OPI40, MET50, PCP10	X
7255F	THC12, COC20, OPI40, AMP50, MET50	
7255ESI	THC12, COC20, OPI40, AMP50, MET50	X

6 Panel

7165ESI	THC12, COC20, OPI40, AMP50, MET50, PCP10	X
7165F	THC12, COC20, OPI40, AMP50, MET50, PCP10	
7264F	THC12, COC20, OPI40, AMP50, MET50, BZO50	
7265F	THC12, COC 20, OPI40, AMP50, MET50, BZO10	
7265ESI	THC12, COC 20, OPI40, AMP50, MET50, BZO10	X

10 Panel



Item #	Configuration & Cut-off	Saliva Saturation Indicator
714105ESI	THC12, COC20, OPI40, AMP50, MET50, PCP10, BZO10, BUP5, OXY20, MTD30	X
714105F	THC12, COC20, OPI40, AMP50, MET50, PCP10, BZO10, BUP5, OXY20, MTD30	
73104F	THC12, COC20, OPI40, AMP50, MET50, PCP10, BZO50, BAR300, BUP10, OXY50	
73105F	THC12, COC20, OPI40, AMP50, MET50, PCP10, BZO10, BAR50, BUP5, OXY50	

11 Panel



74115F-ALC	ALCO, THC12, COC20, OPI40, AMP50, MET50, PCP10, BZO10, BAR50, OXY20, MTD30	
75115ESI	THC12, COC20, OPI40, AMP50, MET50, BZO10, BAR50, BUP5, OXY20, MDMA50, MTD30	X
75115FSI	THC12, COC20, OPI40, AMP50, MET50, BZO10, BAR50, BUP5, OXY20, MDMA50, MTD30	

12 Panel



73125F-ALC	ALCO, THC12, COC20, OPI40, AMP50, MET50, PCP10, BZO10, BAR50, BUP5, OXY20, MTD30	
76125ESI	THC12, COC20, OPI40, AMP50, MET50, PCP10, BZO10, BAR50, BUP5, OXY20, MDMA50, MTD30	X

Please Note: Cut-off levels give approximate detection periods for each substance. The ranges depend on amount and frequency of use, metabolic rate, body mass, age, overall health, drug tolerance, and other additional factors.



STAT Swabs

Item #	Panel	Description
SSB-148	4	THC, COC, MAMP, OPI
SSB-158	5	THC, COC, OPI, MAMP, AMP
SSB-168	6	THC, COC, OPI, MAMP, PCP, AMP
SSB-368	6	THC, COC, OPI, OXY, MAMP, BZO
SSB-168-A	6+ALC	THC, COC, OPI, AMP, MAMP, PCP + ALC
SSB-368-A	6+ALC	THC, COC, OPI, MAMP, BZO, OXY + ALC
SSB-188	8	THC, COC, OPI, MAMP, PCP, AMP, BZO, OXY
SSB-2198-A	9+ALC	THC, COC, OPI, AMP, MAMP, PCP, BZO, MTD, OXY + ALC
SSB-1108	10	THC, COC, OPI, MAMP, PCP, AMP, BZO, OXY, BAR, MTD
SSB-3108	10	THC, COC, OPI, AMP, MAMP, PCP, BZO, MTD, OXY, BUP

Cut-off Levels

Test	Calibrator	Cut-off
Amphetamine (AMP)	D-Amphetamine	50 ng/mL
Methamphetamine (mAMP)	D-Methamphetamine	50 ng/mL
Cocaine (COC)	Benzoylcegonine	20 ng/mL
	Benzoylcegonine	50 ng/mL
Opiate (OPI)	Morphine	40 ng/mL
	Morphine	50 ng/mL
Marijuana (THC)	11-nor- Δ^9 -THC-9-COOH	12 ng/mL
	Δ^9 -THC	25 ng/mL
	Δ^9 -THC	50 ng/mL
Benzodiazepines (BZO)	Oxazepam	50 ng/mL
	Oxazepam	30 ng/mL
Oxycodone (OXY)	Oxycodone	50 ng/mL
Methadone (MTD)	Methadone	75 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Cotinine (COT)	Cotinine	30 ng/mL
Methylenedioxymethamphetamine (MDMA)	(\pm)-3,4-Methylenedioxymethamphetamine	50 ng/mL



MD Saliva Screens

Item #	Configuration	Panel
254	THC, COC, OPI, AMP, MET	5
564	THC, COC, OPI, MET, OXY, BZO	6
274	THC, COC, OPI, MET, OXY, BZO, MDMA	7
874	THC, COC, OPI, MET, AMP, PCP, MDMA	7
6101	THC, COC, OPI, MET, AMP, BUP, MTD, OXY, PCP, BZO	10
3125ALC	THC, COC, OPI, MET, AMP, BUP, BAR, MTD, OXY, PCP, BZO, ALC	12
6125	THC, COC, OPI, MET, AMP, BUP, BAR, MTD, OXY, PCP, BZO, MDMA	12

Test	Calibrator	Cut-off
Amphetamine (AMP)	D-Amphetamine	50 ng/mL
Methamphetamine (mAMP)	D-Methamphetamine	50 ng/mL
Cocaine (COC)	Benzoyllecgonine	20 ng/mL
Opiates (OPI)	Morphine	40 ng/mL
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	12 ng/mL
	Δ^9 -THC	75 ng/mL
Phencyclidine (PCP)	Phencyclidine	10 ng/mL
Benzodiazepines (BZO)	Oxazepam	50 ng/mL
Oxycodone (OXY)	Oxycodone	50 ng/mL
Methadone (MTD)	Methadone	75 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Alcohol (ACL)	Alcohol	> 0.02 % B.A.C